

# Indonesia → Germany

Dog Import Checklist & Quick Guide | TailHarbor | Updated April 2026

## KEY INFORMATION

### COMPLEX ROUTE: Titer test + 90-day wait required

|                            |   |
|----------------------------|---|
| <b>Preparation Time:</b>   | 4-6 months  |
| <b>Microchip:</b>          | ISO 11784/11785 (must be implanted BEFORE rabies vaccination)   |
| <b>Rabies Vaccination:</b> | Required. Minimum 21 days before travel (primary). Keep boosters current.   |
| <b>Titer Test (FAVN):</b>  | REQUIRED. Blood draw 30+ days after vaccination. 90-day wait from draw date.  |
| <b>Documentation:</b>      | Animal Health Certificate (AHC) issued max 10 days before arrival. Single-use.  |
| <b>Breed Restrictions:</b> | Federal ban: Pit Bull Terrier, American Staffordshire Terrier, Staffordshire Bull Terrier, Bull Terrier + crosses. State... |
| <b>Cropping/Docking:</b>   | Dogs with cropped ears or docked tails may face import restrictions in Germany  |
| <b>Insurance:</b>          | Mandatory liability insurance in some regions/states of Germany   |
| <b>Registration:</b>       | Local registration required within 2 weeks of arrival   |
| <b>Dog Tax:</b>            | €90-€150/year typical, up to €600-€1,000/year for dangerous breeds  |

## USEFUL TIPS

- Start the process as early as possible. Delays with lab results or vet appointments can push back your timeline significantly.
- Book your FAVN titer test at an EU-approved lab early. Processing takes 7-14 days, and if results come back below 0.5 IU/ml, you'll need to re-vaccinate and re-test (resetting the 90-day clock).
- Keep your dog's rabies vaccinations current with no gaps. If there's a lapse, the titer test validity may be lost and you'll need to restart the process.
- Research airline pet policies early. Some airlines have limited cargo spots, especially in summer. Temperature embargoes apply above 30C / below -5C at any transit airport.
- Start crate training weeks before travel. A calm, crate-confident dog will have a much less stressful journey.

## PRE-TRAVEL CHECKLIST

### EARLY PREPARATION

- Confirm your dog's breed is permitted in Germany
- Research airline pet policies and book cargo/cabin spot
- Purchase IATA-approved travel crate (correct size for your dog)
- Begin crate training

### MEDICAL REQUIREMENTS

- Microchip implanted (ISO 11784/11785) -- BEFORE vaccination
- Rabies vaccination administered (min. 21 days before travel)
- Blood draw for FAVN titer test (30+ days after vaccination)
- Titer test result received (minimum 0.5 IU/ml)
- 90-day waiting period completed (from blood draw date)

### DOCUMENTATION

- Animal Health Certificate obtained (max 10 days before arrival)
- All vaccination records compiled
- Titer test certificate included
- Copies of all documents made (2 sets)
- Airline booking confirmation printed

### TRAVEL DAY

- All original documents in carry-on bag
- Crate labeled with your name, phone, destination, and 'Live Animal'
- Water dish filled and attached inside crate
- Arrive at airport early (3+ hours for international with pet)

### AFTER ARRIVAL IN GERMANY

- Present documents at Travellers' Point of Entry / Border Control
- Register dog with local municipality within 2 weeks
- Register with local veterinarian
- Check regional insurance requirements and obtain liability insurance if required